



Analytical[®]
Technologies Limited



Products | Applications | Software | Accessories | Consumables | Services

Analytical Technologies Limited
www.analyticalgroup.net

SPECIFICATION

CCD camera : Professional Colour digital CCD Camera
Pixels : 1.4 million
AD converter : 12 bit
Camera Lens : 8mm Megapixel lens
Illumination : Long-life upper incident
and bottom transmitted LED arrays
Petri dish Diameter : Standard 50mm~110mm
Testing method : Fully automatic with manual option
Counting time : Over 300 colonies per
second, 1 plate per second
Detectable Limit : Colony 0.1mm
Zone 0.05mm
Data interface : USB2.0
Software : Users-friendly software with colony
counting & zone measuring function



Analytical®
Technologies Limited

Analytical Group Companies
Analytical Technologies Limited
Analytical Bio-Technologies
Analytical Bio-Med
HPLC Solutions

Branch Office / Distributor

Corporate & Regd. Office :

Analytical House, # E67 & E68. Ravi Park, Opp., Sainath Park, Vasna Road, Baroda, Gujarat - 390 015. INDIA
Tele : 0265-2253620 / 2252839 Fax : 0265-2254395 E-mail : info@hplctechnologies.com / info@ais-india.com
Web : www.ais-india.com / www.analyticalgroup.net