



Pensky Martin Flash Point Apparatus (Model ATL-FPP)



High quality components are used in the **Pensky Martin Flash Point Apparatus** that we offer to clients. The Apparatus is designed with accuracy which ensures smooth and flawless performance. We are reputed as the foremost company that brings reliable Pensky Martin Flash Point Apparatus to the clients.

Application

Determining close cup flash points of oils, and other viscous material and the suspension of solids

WITH DIGITAL TEMP. CONTROLLER

Pensky Marten Flash Point App.

Electrical Heating with E.R. fitted with motorised FHP elect. Stirrer with stand. gas Test jet/ Oil Test jet With Digital Temperature Indicator cum Controller

»» About Analytical Technologies

Analytical Technologies is a Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO : 9001 Certified Company engaged in Designing, Manufacturing, Marketing & Providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, & Clinical Diagnostics Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers.

»» Our Other Products And Technologies



NOVA 2030 Plus
Automatic Bio-Chemistry
Analyzer



HEMA 1062 Plus
Fully-automatic
Hematology Analyzer



NOVA 2200 Plus
Auto
Chemistry Analyzer



URINOVA 2800
Urine Analyzer



Fully Automated
CLIA



Micro Plate
Reader / Washer



MRI CT Scanner



NOVALYTE 2200
Electrolyte Analyzer



NOVA Basic Semi-Auto
Chemistry Analyzer



Kidney Strain
Analyzer



Amino Acid
Analyzer



PCR/Gradient PCR/RT-PCR



5 part
Hematology Analyzer



Blood Gas Analyzer



Portable POC
Blood gas Analyzer



Random access analyzer
for immunoassay,
proteins & clinical chemistry



Revolutionary automation
in neonatal screening



SPECTRO 2092



Semen Analyzer



Water Purification
System



In Chemistry
Reagents



Urine Strips



Rapid Diagnostic kits



Hematology
Reagents



Analytical[®]
Technologies Limited

Analytical Group Companies
Analytical Technologies Limited
Analytical Bio-Technologies
Analytical Bio-Med
HPLC Solutions

Corporate & Regd. Office :

Analytical House, # E67 & E68, Ravi Park, Opp., Sainath Park, Vasna Road, Baroda, Gujarat - 390 015, INDIA

Tele : 0265-2253620 / 2252839 Fax : 0265-2254395

E-mail : info@hplctechnologies.com / info@ais-india.com Web : www.ais-india.com / www.analyticalgroup.net