



## Analytical Technologies Fetoscope



### ATL 3004 FETO

#### Specification:-

Handheld light weight pocket - size device based on Ultrasonic Doppler to measure foetal heart rate (FHR) and listen to the sound of the Foetal heart beat.

#### A. Specifications

1. 2MHz or 3MHz waterproof probe
2. FHR Range : 50 – 200 BPM , resolution : 1 BPM
3. LCD / Digital Display of the heart beat rate
4. Power source : Rechargeable AA Batteries
5. Built-in loudspeaker and adjustable volume listening to the baby sound through headphones
6. Headphone output
7. Headphone set
8. Battery charger
9. Carrying Case