



INCUBATOR BACTERIOLOGICAL



Products | Applications | Software | Accessories | Consumables | Services

INCUBATOR BACTERIOLOGICAL

Standard double wall fabrication. Inner chamber made out of richly anodized aluminium or highly polished stainless steel sheet and outer made out of thick mild steel sheet finished in white stoving enamel/powder coated paint with mat finished colour combinations. Double wall door with double glass window for observation in the chamber. Temperature range 5° C above ambient to 80° C ± 1° C controlled by a thermostat. The equipment is workable on 220 VAc 50 Hz single phase.

Chamber size in mm/Inches	No. of	Capacity
H x D x W	Shelves	(Approx.)
600 x 450 x 450 (24"x 18"x18")	2	125 ltrs.



Analytical Group Companies
Analytical Technologies Limited
Analytical Bio-Technologies
Analytical Bio- Med
HPLC Solutions

Corporate & Regd. Office :

Analytical House, # E67 & E68, Ravi Park, Opp., Sainath Park, Vasna Road, Baroda, Gujarat - 390 015. INDIA

Tele : 0265-2253620 / 2252839 Fax : 0265-2254395

E-mail : info@hplctechnologies.com / info@ais-india.com Web : www.ais-india.com / www.analyticalgroup.net