



**Analytical**<sup>®</sup>  
Technologies Limited

## KJELDAHL DISTILLATION AND DIGESTION UNIT



### SPECIFICATION

This unique combined unit is supplied with two sets of Mantle type heating platforms, one for distillation and one for digestion. The unit is also provided with condenser rack, lead fume duct covered by S.S. frame and pair of hangers. Complete set is supplied on a M.S. Tubular stand. Both units are supplied with individual energy regulators. The Unit is of 3/6/test Mantle/ Heater Type and work on 220/ 230 volts A.C. supply.



**Analytical**<sup>®</sup>  
Technologies Limited

Analytical Group Companies  
Analytical Technologies Limited  
Analytical Bio-Technologies  
Analytical Bio- Med  
HPLC Solutions

Branch Office / Distributor

Corporate & Regd. Office :

Analytical House, # E67 & E68. Ravi Park, Opp., Sainath Park, Vasna Road, Baroda, Gujarat - 390 015. INDIA  
Tele : 0265-2253620 / 2252839 Fax : 0265-2254395 E-mail : info@hplctechnologies.com / info@ais-india.com  
Web : www.ais-india.com / www.analyticalgroup.net