



# AMIA-S

## Standard Milk Analyzer

Model	AMIA-S		
Measuring time	Measuring parameters	Standard complete	Additional options
90 Sec	FAT – 0.01% - 25% SNF – 3% - 15% Density – 1015 – 1040kg/m <sup>3</sup> Proteins – 2% - 7% Lactose – 0.01% - 6% Water content – 0% - 70% Temp. of milk – 1°C– 40°C Freezing point - -0.400 – 0.700°C Solids – 0.4 – 1.5%	1. Sample visualization on LCD display(green) 2. Possibility to measure 1 sample – 90 seconds (option 60 seconds) 3. Adapter 220-110 V 4. Spare silicone tubes 5. Syringe 6. Plastic sample holder 7. Calibrations Cow, Sheep, and UHT milk 8. Operation Manual 9. Cardboard transport case 10. CD – Service Pack	Ph option
60 Sec			Conductivity
	High Fat 45% (High Density)		

## ▶▶ About Analytical Technologies

Analytical Technologies is a Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO : 9001 Certified Company engaged in Designing, Manufacturing, Marketing & Providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, & Clinical Diagnostics Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers.

## ▶▶ Our Other Products And Technologies



NOVA 2020 Plus  
Automatic Bio-Chemistry  
Analyzer



HEMA 2062 Plus  
Fully-automatic  
Hematology Analyzer



NOVA 2200 Plus  
Auto  
Chemistry Analyzer




URINOVA 2800  
Urine Analyzer



Fully Automated  
CLIA



Micro Plate  
Reader / Washer



MRI CT Scanner



NOVALYTE 2200  
Electrolyte Analyzer



NOVA Basic Semi - Auto  
Chemistry Analyzer



Kidney Stone  
Analyzer



Amino Acid  
Analyzer



PCR/Gradient PCR/RT PCR



5 part  
Hematology Analyzer



Blood Gas Analyzer



Portable POC  
Blood gas Analyser



Random access analyzer  
for immunoassay,  
proteins & clinical chemistry



Revolutionary automation  
in neonatal screening



SPECTRO 2092



Semen Analyzer



Water Purification  
System



Bio Chemistry  
Reagents



Urine Strips



Rapid Diagnostic kits



Hematology  
Reagents



**Analytical**®  
Technologies Limited

Analytical Group Companies  
Analytical Technologies Limited  
Analytical Bio-Technologies  
Analytical Bio- Med  
HPLC Solutions

Corporate & Regd. Office :

Analytical House, # E67 & E68. Ravi Park, Opp., Sainath Park, Vasna Road, Baroda, Gujarat - 390 015. INDIA

Tele : 0265-2253620 / 2252839 Fax : 0265-2254395

E-mail : [info@hplctechnologies.com](mailto:info@hplctechnologies.com) / [info@ais-india.com](mailto:info@ais-india.com) Web : [www.ais-india.com](http://www.ais-india.com) / [www.analyticalgroup.net](http://www.analyticalgroup.net)