



Rotary Vacuum Evaporator

4 Litre Rotary Evaporator

Specification :

Evaporating Flask	:- 4000ml
Receiving Flask	:- 1000ml
Cooling Area	:- 0.13 M2
Rotary Speed	:- 20-300rpm
Power of Motor	:- 0.40 kw
Bath Power	:- 1.5KW
Bath Lift	:- 1000mm - Manual
Sealing System	:- Rubber vacuum seal + PTFE Sleeve
Temperature Control Range	:- 0-180°C (Intelligent Digital Control)
Temperature Wave Range	:- +/- 0.2°C
Dimension	:- 620 x 330 x 820 mm
Power	:- 220V/50Hz

▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

▶▶▶ Reach us @



HPLC Solutions		MultipleLabs		Analytical Bio-Med		Analytical Distributors		Analytical Foundation (Trust)	
Corporate & Regd. Office:		T: +91 265 2253620		E: info@hplctechnologies.com		W: www.ais-india.com		Sales & Support Offices:	
Analytical House, # E67 & E68,		+91 265 2252839		info@multiplelabs.com		www.analyticalgroup.net		across the country :	
Ravi Park, Vasna Road, Baroda,		+91 265 2252370		info@analyticalgroup.net		www.hplctechnologies.com		Distributors & Channel	
Gujarat 390 015. INDIA		F: +91 265 2254395		info@analyticalbiomed.com		www.multiplelabs.com		partners World Wide	

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.