



Video Bronchoscope (Paediatric) BRONCHO 3700 Series

Technical Specification:

Diagnostic Laryngo /Pharyngoscope with Fiber Optic Cable & Light Source One Biopsy Cable Should Be Provided With The Set

- Field of View- 120 degree
- Depth of field- 3-100mm;
- Distal end Outer Dia: 3.8 mm or less,
- Insertion Tube, Outer Dia: 3.8mmor less;

- Bending section : Angulation Up 180 degree, Down 130 degree; Working Length 600mm or more;
- Minimum visible distance 3mm. Inner diameter of instrument channel- 1.2mm or less.
- Include the following accessories: Biopsy forceps- 1, Cytology brush- 1 Pk; Cannula standard type- 1; Mouth piece- 2; Semi-disposable biopsy valve- 30 pcs; Semi disposable suction valve- 10 pcs; Cleaning and maintenance kit- 1

We are counted amongst the apex organizations involved in manufacturing, importing and supplying qualitative **Fiber Optic Bronchoscope**. In our advanced manufacturing unit, this product is manufactured with ultra-grade raw material implementing sophisticated technologies to ensure its longer service life at the clients' end. We stringently tested it on different parameters so as to make it comply with the international norms.

Features:

Customized length
 High temperature resistance
 Excellent performance
 Cost effective



Analytical®

Technologies Limited

HPLC Solutions Analytical Bio- Med MultipleLabs

Analytical House, # E67 & E68, Ravi Park, Opp. Sainath Park, Vasna Road, Baroda - 390 015. Gujarat, INDIA.
 Tele.: +91 265 225 3620, 225 2839 Fax: +91 265 225 4395 Web: www.analyticalgroup.net / www.ais-india.com
 Email: info@ais-india.com / info@hplctechnologies.com / info@multiplelabs.com