



**Analytical**<sup>®</sup>  
Technologies Limited



## **Digital Trolley Ultrasound Scanner**

Products | Applications | Software | Accessories | Consumables | Services

**Analytical Technologies Limited**  
[www.analyticalgroup.net](http://www.analyticalgroup.net)

## Trolley Digital Ultrasound Scanner

Digital trolley ultrasound, economic type, 8-segment TGC Main unit, 3.5MHz convex probe(80 element), 14" SVGA monitor, peripheral ports of video and SVGA

### Main features

- Full digital with latest ultrasonic imaging techniques
- Back-light, silicon-gel alphanumeric keyboard
- Applicable for abdomen, GYN, OB, cardiac, urology and small parts
- 5 working modes: B, B/B, B/M, M and 4B
- 12" flicker free SVGA monitor, backlit silicon-gel keyboard
- 2 probe sockets
- 8 segment TGC
- Permanent storage of 32 frames of ultrasound images
- Peripheral ports of USB, RS232, video, and grounding

### Standard configurations

- CUS-9618C main unit
- 12" flash-free monitor
- 3.5MHz convex probe
- Coupling gel
- USB drive disk and cable

**Products | Applications | Software | Accessories | Consumables | Services**

**Analytical Technologies Limited**  
[www.analyticalgroup.net](http://www.analyticalgroup.net)